GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™

A Rapid Membrane Enzyme Immunoassay for the Simultaneous Qualitative Detection and Differentiation of *Giardia* Cyst Antigen and *Cryptosporidium* Oocyst Antigen in Human Fecal Specimens

Catalog No. T30407 (25 Tests)

IVD *In Vitro* Diagnostic Medical Device

For Canadian Users: For Laboratory Use Only

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GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™

INTENDED USE

The *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test is a rapid membrane enzyme immunoassay for the simultaneous qualitative detection and differentiation of *Giardia* cyst antigen and *Cryptosporidium* oocyst antigen in a single test device. It is intended for use with human fecal specimens from patients with gastrointestinal symptoms to aid in the diagnosis of *Giardia* and/or *Cryptosporidium* gastrointestinal infection. The test results should be considered in conjunction with the patient history.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

EXPLANATION

Giardia spp. is a binucleated flagellated protozoan parasite which exists in two forms: a noninfectious, pear-shaped trophozoite (9 to 20 µm) inhabiting the small intestine, and the highly infectious cyst form which is elliptical in shape and ranges in size from 8 to 12 µm (1). Survival outside its host varies greatly between the two forms: the trophozoite which is extremely labile, lasts only a matter of hours outside the body, while the cyst form may survive for several days in external environments (1). Travelers have been found to contract giardiasis from endemic areas, specifically contaminated water (2-5). Transmission also occurs by direct contact especially by asymptomatic carriers and by food contamination (6-7). High-risk individuals include young children, immunocompromised patients and those without previous exposure (8). More recently, giardiasis has become a common sexually transmitted disease (9). Animal fecal contamination especially of water is another route of transmission in humans (4, 10, 11). Clinical manifestations of giardiasis range from asymptomatic carriage with the passing of cysts to chronic debilitating diarrhea, weight loss and malabsorption (8,12,13).

Cryptosporidium spp. is a protozoan parasite of vertebrates previously thought to cause diarrhea only in animals (14). In 1976, the first human infection was reported (15). Subsequently it has been found to be associated with diarrheal illness world-wide and is a frequent cause of traveler's diarrhea (14,16). The disease is transmitted by the thick-walled oocyst form, 2-6µm in diameter, which is remarkably resistant to common disinfectants and routine chlorination of drinking water. Person-to-person transmission, especially among children, is common (17). *Cryptosporidium* has little to no host specificity and animals such as rodents, cattle, and domestic pets serve as a reservoir for zoonotic transmission to humans (14,18). This occurs either by direct contact or by contamination of water supplies with fecal matter (14,19-21). Cryptosporidiosis is a serious opportunistic infection in patients with acquired immunodeficiency syndrome (AIDS) and is potentially a sexually transmitted disease (19,22). Clinical manifestations of cryptosporidiosis include cholera-like diarrhea, abdominal pain, nausea, vomiting, and weight loss. In normal persons the infection is usually self-limiting and short term. In AIDS and immunosuppressed patients, cryptosporidiosis can result in prolonged and life-threatening illness due to excessive fluid loss. In these patients the infection may also spread to the respiratory and biliary tracts (19).

PRINCIPLE OF THE TEST

The *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test uses monoclonal and polyclonal antibodies to cell-surface antigens of the organisms. The device contains a *Reaction Window* with three vertical lines of immobilized antibodies. The *Giardia* test line ("Giar") contains mouse monoclonal antibodies against *Giardia*. The *Crypto* test line ("Cryp") contains mouse monoclonal antibodies against *Cryptosporidium*. The control line ("C") is a dotted line that contains anti-horseradish peroxidase (HRP) antibodies. The *Conjugate* consists of polyclonal antibodies coupled to horseradish peroxidase. To perform the test, the sample is added to a tube containing a mixture of *Diluent* and *Conjugate*. The diluted sample-conjugate mixture is added to the *Sample Well* and the device is allowed to incubate at room temperature for 15 minutes. During the incubation, cyst and/or oocyst

antigens in the sample bind the antibody-peroxidase conjugates. The antigen-antibody-conjugate complexes migrate through a filter pad to a membrane where they are captured by the immobilized *Giardia* and/or *Cryptosporidium*-specific antibodies in the test lines. The *Reaction Window* is subsequently washed with *Wash Buffer*, followed by the addition of *Substrate*. After a 10 minute incubation period, the reaction is examined visually for the appearance of a vertical blue line on either side of the *Reaction Window*. A blue line indicates a positive test. A positive “control” reaction, indicated by a vertical dotted blue line under the “C” portion of the *Reaction Window*, confirms that the test is working properly and the results are valid.

MATERIALS PROVIDED

MEM DEV	Membrane Devices – each pouch contains 1 device
DIL SPE	Diluent (20 mL) – Buffered protein solution with graduated dropper assembly*
WASH REAG	Wash Buffer (12 mL) – Buffered solution with graduated dropper assembly*
SUBS REAG	Substrate (3.5 mL) – Solution containing tetramethylbenzidine
CONJ ENZ	Conjugate (2.0 mL) – Polyclonal antibody to a cell-surface antigen of <i>Giardia</i> coupled to horseradish peroxidase and a polyclonal antibody to a cell-surface antigen of <i>Cryptosporidium</i> coupled to horseradish peroxidase in a buffered protein solution*
CONTROL +	Positive Control (2.5 mL) – Inactivated <i>Giardia</i> and <i>Cryptosporidium</i> antigens in a buffered protein solution*

Disposable plastic pipettes – graduated at 25 µL, 100 µL, 200 µL, 300 µL, 400 µL, and 500 µL

*contains 0.05% ProClin® 300

Signal Word: Warning

H317: May cause an allergic skin reaction

P261, P272, P280, P302, P352, P333, P313, P321, P362, P363, P364, P501

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

<i>Small test tubes (e.g., plastic Eppendorf tubes)</i>	<i>Applicator sticks</i>
<i>Timer</i>	<i>Vortex mixer</i>
<i>Disposable gloves for handling fecal samples</i>	<i>Pipettor and tips</i>

SHELF LIFE AND STORAGE

The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. The kit should be stored between 2° and 8°C.

PRECAUTIONS

1. Rx - Prescription Only
2. For *in vitro* diagnostic use. For professional use only.
3. Each component in the kit should be inspected for any signs of leakage. Upon arrival, inspect the kit to ensure that components are not frozen or warm to the touch due to improper shipping conditions.
4. The *Substrate* reagent should be colorless. If the *Substrate* reagent changes to a dark blue/violet color, discard and call technical services for replacement.
5. Reagents from different kits should not be mixed. Do not use a kit past the expiration date.
6. Caps, tips, and dropper assemblies are color-coded; do NOT mix or interchange!
7. Bring all components to ROOM TEMPERATURE BEFORE USE!
8. Do not freeze the reagents. The kit should be stored between 2°C and 8°C.
9. The pouch containing the *Membrane Device* should be at room temperature before opening. Keep the *Membrane Devices* dry before use.
10. Hold reagent bottles vertically to dispense reagents to ensure consistent drop size and correct volume.
11. Microbial contamination of reagents may decrease the accuracy of the assay. Avoid microbial contamination of reagents by using sterile disposable pipettes if removing aliquots from reagent bottles.
12. *Membrane Devices* cannot be reused.

13. The test has been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test. Do not deviate from the specified procedure.
14. Be attentive to the total assay time when testing more than one fecal specimen. Add *Diluent* first, and then add the *Conjugate* to each tube of *Diluent*. Then add specimen to the tube of *Diluent/Conjugate*. Thoroughly mix all of the diluted specimens, and transfer to the *Membrane Device*. The 15-minute incubation step begins after the last diluted sample-conjugate mixture has been transferred to the *Membrane Device*.
15. Do not concentrate specimens before testing.
16. Use fresh fecal specimens within 72 hours of collection to obtain optimal results. Specimens that are frozen may lose activity due to freezing and thawing. If using frozen specimens, thaw at room temperature.
17. This test has been shown to be compatible with specimens preserved in 10% formalin, sodium acetate formalin, and transport media such as Cary Blair or C&S when following the "Sample Preparation" guidelines listed in this Package Insert. However, the test has not been proven compatible with other preservatives and transport media.
18. Specimens and *Membrane Devices* should be handled and disposed of as potential biohazards after use. Wear disposable gloves when doing the test.
19. Fecal specimens may contain potentially infectious agents and should be handled at "Biosafety Level 2" as recommended in the CDC/NIH Manual, "Biosafety in Microbiological and Biomedical Laboratories."
20. The *Conjugate*, *Diluent*, *Positive Control*, and *Wash Buffer* reagents contain 0.05% ProClin® 300 as a preservative. Although the concentration is low, ProClin® 300 is known to be harmful (skin sensitization may occur). If skin sensitization/irritation or rash occurs, get medical advice/attention. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
21. Follow your national, regional, and local ordinances accordingly for waste disposal regulations.

COLLECTION AND HANDLING OF FRESH FECAL SPECIMENS

Specimens collected for routine ova and parasite examination can be used with the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test. Fecal specimens should be collected in clean, leak-proof containers.


1. Standard collection and handling procedures used in-house for fecal specimens are appropriate.
2. Fresh, untreated specimens should be stored between 2° and 8°C. Test fresh specimens that are less than 72 hours old, whenever possible. Store fresh specimens frozen ($\leq -10^{\circ}\text{C}$) for up to 90 days if the test cannot be performed within 72 hours of collection, but note that freezing and thawing of the specimen may result in loss of activity due to degradation of the antigen. Avoid multiple freeze-thaw cycles. If using frozen specimens, thaw at room temperature.
3. Fecal specimens that have been concentrated or treated with PVA fixatives are not suitable for use with this test.
4. Make sure that specimens are thoroughly mixed PRIOR to performing the assay.
5. Storing fecal specimens in the *Diluent* is NOT recommended.

SAMPLE PREPARATION

1. Bring all reagents and the required number of devices to room temperature before use.
2. Set up and label one small test tube for each specimen and optional external controls as necessary.
3. Using the gray graduated dropper assembly, add 500 μL (2nd graduation from the tip) *Diluent* to each tube for fecal specimens. For specimens in fixative or transport media such as 10% formalin, add 400 μL of *Diluent* to the tube. For external controls add 400 μL of *Diluent* to each tube.

4. Add one drop of *Conjugate* (red capped bottle) to each tube.
5. Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample – the pipettes have raised graduations at 25 μL , 100 μL , 200 μL , 300 μL , 400 μL , and 500 μL .

Graduated Transfer Pipette:

6. Mix all specimens thoroughly regardless of consistency – it is essential that the specimens be evenly suspended before transferring.

Liquid/Semi-solid specimens – pipette 25 μL of specimen with a transfer pipette (graduated at 25 μL , 100 μL , 200 μL , 300 μL , 400 μL , and 500 μL) and dispense into the *Diluent/Conjugate* mixture. Use the same transfer pipette to mix the diluted specimen.

Formed/Solid specimens – Care must be taken to add the correct amount of formed feces to the sample mixture. Mix the specimen thoroughly using a wooden applicator stick and transfer a small portion (approximately 2 mm diameter, the equivalent of 25 μL) of the specimen into the *Diluent/Conjugate* mixture. Emulsify the specimen using the applicator stick.

Fecal specimens in fixative or transport media – pipette 100 μL (2 drops from transfer pipette) of sample into the *Diluent/Conjugate* mixture.

7. **Optional External Control Samples:**

External Positive Control – add four drops of *Positive Control* (gray-capped bottle) to the appropriate test tube containing *Diluent*.


External Negative Control - add 100 μL *Diluent* to the appropriate test tube containing *Diluent*.

NOTE: Transferring too little specimen, or failure to mix and completely suspend the specimen in the *Diluent/Conjugate* mixture, may result in a false-negative test result. The addition of too much fecal specimen may cause invalid results due to restricted sample flow.

TEST PROCEDURE

1. Obtain one *Membrane Device* per specimen, and one *Membrane Device* per optional external positive or negative control as necessary. The foil bags containing the devices should be brought to room temperature before opening. Use the device immediately after opening. Label each device appropriately and orient it on a flat surface so the “QUIK CHEK” print is located at the bottom of the device, and the small *Sample Well* is located in the top right corner of the device.

Membrane Device

2. Close each tube of diluted specimen and mix thoroughly. Proper mixing can be achieved by vortexing or inverting the tube. Once a patient sample or *Positive Control* has been diluted in the *Diluent/Conjugate* mixture, it may be incubated at 2 - 8°C for any period of time up to 24 hours prior to addition to the *Membrane Device*.
3. Using a new transfer pipette, transfer 500 μL of the diluted sample-conjugate mixture into the **Sample Well** (smaller hole in the top right corner of the device) of a *Membrane Device*, making certain to expel the liquid sample onto the wicking pad inside of the *Membrane Device*. When loading the sample into the sample well, make sure that the tip of the transfer pipette is angled towards the *Reaction Window* (larger hole in the middle of the device).

4. Incubate the device at room temperature for 15 minutes – the sample will wick through the device and a wet area will spread across the *Reaction Window*.

NOTE FOR SAMPLES THAT FAIL TO MIGRATE:

Occasionally, a diluted fecal specimen clogs the membrane and the Reaction Window does not wet properly. This is ordinarily due to the addition of too much fecal specimen to the sample Diluent. If the diluted fecal specimen fails to migrate properly within 5 minutes of adding the sample to the Sample Well (i.e. the membrane in the Reaction Window does not appear to be completely wet), then add 100 µL (4 drops) of Diluent to the Sample Well and wait an additional 5 minutes (for a total of 20 minutes). If the specimen still fails to migrate, retest the specimen.

5. After the incubation, add 300 µL of *Wash Buffer* to the **Reaction Window** using the graduated white dropper assembly (or equivalent). Allow the *Wash Buffer* to flow through the *Reaction Window* membrane and be absorbed completely.
6. Add 2 drops of *Substrate* (white-capped bottle) to the **Reaction Window**. Read and record results visually after 10 minutes.

INTERPRETATION OF RESULTS

1. Interpretation of the test is most reliable when the device is read immediately at the end of the reaction period. Read the device at a normal working distance in a well-lit area. View with a line of vision directly over the device. A positive result may be interpreted at any time between the addition of *Substrate* and the 10-minute read time. However, a test cannot be interpreted as negative or invalid until the 10 minutes following the addition of the *Substrate* has been completed.
2. Observe the device for the appearance of blue dots in the middle of the *Reaction Window* representing the internal positive control. Observe device for the appearance of blue lines on either side of the control dots (*Giardia* on the right, *Cryptosporidium* on the left). Lines may appear faint to dark in intensity.
3. **Positive *Giardia* Result (Fig. 1b, Fig. 1c):** The control dots and the *Giardia* line (on the right side of the *Reaction Window*) are visible. The lines may appear faint to dark in intensity - any blue line on the right side of control dots is considered positive. Do not interpret membrane discoloration as a positive result. A positive result indicates the presence of *Giardia* antigen and a properly reactive positive control line.
4. **Positive *Cryptosporidium* Result (Fig. 1b, Fig. 1d):** The control dots and the *Cryptosporidium* line (on the left side of the *Reaction Window*) are visible. The lines may appear faint to dark in intensity - any blue line on the left side of control dots is considered positive. Do not interpret membrane discoloration as a positive result. A positive result indicates the presence of *Cryptosporidium* antigen and a properly reactive positive control line.
5. **Negative Result (Fig. 1a):** Only the control dots in the middle of the *Reaction Window* are visible. No lines are visible. A negative result indicates that *Giardia* and *Cryptosporidium* antigen are either not present in the sample, or are below the detectable limits of this test.
6. **Invalid Result (Fig. 1e-1h):** The test is invalid if the control dots are not present at the completion of the reaction period (Figures 1e-1h).

FIGURE 1: GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ INTERPRETATION OF RESULTS


Figure 1a
Negative Result

Figure 1b
Positive *Giardia* and
Cryptosporidium Result

Figure 1c
Positive *Giardia*
Result

Figure 1d
Positive *Cryptosporidium*
Result


Figure 1e
Invalid Result


Figure 1f
Invalid Result


Figure 1g
Invalid Result


Figure 1h
Invalid Result

QUALITY CONTROL

Internal: A dotted blue line must be visible in the middle of the *Reaction Window*, below the “C” on every *Membrane Device* that is tested. The appearance of the blue control dots confirms that the sample and reagents were added correctly, that the reagents were active at the time of performing the assay, and that the sample migrated properly through the *Membrane Device*. A clear background in the result area is considered an internal negative control. If the test has been performed correctly and reagents are working properly, the background will be white to give a discernable result.

External: The reactivity of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* kit should be verified upon receipt using the *Positive Control* and negative control (*Diluent*). The *Positive Control* is supplied with the kit (gray-capped bottle). The *Positive Control* confirms the reactivity of the other reagents associated with the assay. *Diluent* is used for the negative control. Additional tests can be performed with the controls to meet the requirements of local, state, and/or federal regulations and/or accrediting organizations.

LIMITATIONS

1. The *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test is used to detect *Giardia* and/or *Cryptosporidium* in fecal specimens. The test confirms the presence of antigen in feces and this information should be taken under consideration by the physician in light of the clinical history and physical examination of the patient.
2. Fecal specimens are extremely complex. Optimal results with the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test are obtained with specimens that are less than 72 hours old. Most undiluted fresh specimens can be stored between 2°C and 8°C for 72 hours before degradation of the antigen is noted. If specimens are not assayed within this time period, they may be frozen for up to 90 days and thawed. However, repeated freezing and thawing may result in loss in the immunoreactivity of antigen.
3. Some specimens may give weak reactions. This may be due to a number of factors such as the presence of low levels of antigen, the presence of binding substances, or inactivating enzymes in the feces. *Under these conditions, a fresh specimen should be tested.*
4. The *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test is qualitative. The intensity of the color should not be interpreted quantitatively.
5. The use of concentrated fecal specimens has not been evaluated.

EXPECTED VALUES

Normal healthy individuals should not be infected with *Giardia* or *Cryptosporidium* and should test negative in the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test. A positive test result in the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test indicates that the person is shedding detectable amounts of *Giardia* and/or *Cryptosporidium* antigen. The incidence of *Giardia* and *Cryptosporidium* infection varies significantly between populations and geographic regions. Children in daycare settings have exhibited higher rates of infection with *Giardia* than the normal population (21). In addition, homosexual men have shown higher rates of infection (22). In general, laboratory-confirmed incidence of cryptosporidiosis in developed countries ranges from 1 to 2% overall with a higher incidence in children (23).

Prospective studies to evaluate the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test were conducted at an Asian hospital specializing in diarrheal diseases. The studies included 129 freshly collected fecal specimens. A positive test result in the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test was obtained in 22.5% for *Giardia* and in 5.4% for *Cryptosporidium* antigen.

PERFORMANCE CHARACTERISTICS

The performance of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test was evaluated at 3 geographically diverse sites. At Site #1 and Site #3 the performance of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test was compared to Microscopy (IFA) and included 220 fresh, 140 frozen, 216 preserved-formalin, and 215 preserved-SAF. At Site #1 and Site #2 the performance of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test was compared to two commercially available ELISAs (predicate devices for *Giardia* and *Cryptosporidium*) and included 349 fresh, 322 frozen, 36 preserved-formalin, and 142 preserved-SAF.

Performance as compared to Microscopy (IFA) - Combined Results for Study Sites #1 and #3:

For *Giardia* spp.

The following table shows a summary of the clinical performance of the *Giardia* portion of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test. The results show that the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test exhibited a sensitivity of 98.9%, a specificity of 100%, and an overall correlation of 99.7% with Microscopy – IFA (considered the gold standard).

Combined Results - Clinical Performance Comparing the *Giardia* Line of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* Test to Microscopy – IFA

N = 791	Microscopy - IFA <i>Giardia</i> positive	Microscopy - IFA <i>Giardia</i> negative
<i>GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™</i> - <i>Giardia</i> Line Positive	181	0
<i>GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™</i> - <i>Giardia</i> Line Negative	2	608

		95% Confidence Limits
Sensitivity	98.9%	95.7 – 99.8%
Specificity	100%	99.2 – 100%
Correlation	99.7%	99.7 – 99.7%

For *Cryptosporidium* spp.

The following table shows a summary of the clinical performance of the *Cryptosporidium* portion of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test. The results show that the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test exhibited a sensitivity of 100%, a specificity of 99.8%, and an overall correlation of 99.9% with Microscopy – IFA (considered the gold standard).

Combined Results - Clinical Performance Comparing the *Cryptosporidium* Line of the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* Test to Microscopy – IFA

N = 791	Microscopy - IFA <i>Crypto.</i> positive	Microscopy - IFA <i>Crypto.</i> negative
<i>GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™</i> - <i>Cryptosporidium</i> Line Positive	140	1
<i>GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™</i> - <i>Cryptosporidium</i> Line Negative	0	650

		95% Confidence Limits
Sensitivity	100%	96.7 – 100%
Specificity	99.8%	99.0 – 100%
Correlation	99.9%	100 – 100%

Performance as Compared to Two Predicate Devices (commercially available devices for *Giardia* and *Cryptosporidium*) - Combined Results for Study Site #1 and Site #2:

The combined results of our performance evaluations at Study Site #1 and #2 as compared to two commercially available ELISAs (predicate devices for *Giardia* and *Cryptosporidium*) exhibited a 99.1% agreement for *Giardia* positive specimens, a 99.7% agreement for *Giardia* negative specimens, with an overall agreement of 99.5%. The test exhibited a 99.2% agreement for *Cryptosporidium* positive specimens, 99.6% agreement for *Cryptosporidium* negative specimens, and an overall agreement of 99.5%.

***GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test *Giardia* Line versus the Commercial ELISA for *Giardia* Detection**

N = 849	Commercial ELISA - <i>Giardia</i> positive	Commercial ELISA - <i>Giardia</i> negative
<i>GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™</i> - <i>Giardia</i> Line Positive	213	2
<i>GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™</i> - <i>Giardia</i> Line Negative	2	632

		95% Confidence Limits
Sensitivity	99.1%	96.3 – 99.8%
Specificity	99.7%	98.7 – 99.9%
Correlation	99.5%	99.5 – 99.5%

GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ test *Cryptosporidium* Line versus the Commercial ELISA for *Cryptosporidium* Detection

N = 849	Commercial ELISA - <i>Cryptosporidium</i> positive	Commercial ELISA - <i>Cryptosporidium</i> negative
GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ - <i>Cryptosporidium</i> Line Positive	130	3
GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ - <i>Cryptosporidium</i> Line Negative	1	715

		95% Confidence Limits
Sensitivity	99.2%	95.2 – 100%
Specificity	99.6%	98.7 – 99.9%
Correlation	99.5%	99.5 – 99.5%

ANALYTICAL SENSITIVITY

The analytical sensitivity of the device was determined by spiking purified *Giardia* cysts or *Cryptosporidium* oocysts quantified by immunofluorescent antibody microscopy (IFA) into negative human fecal specimens. The concentration of *Giardia* cysts and *Cryptosporidium* oocysts in fecal matrix where specimens were positive by the GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ test 95% of the time were utilized to describe the assay limit-of-detection (LOD). Test results determined the LOD for the assay to be 6000 cysts/mL of feces for *Giardia* (equivalent to 133 cysts detected per test) and 6000 oocysts/mL feces for *Cryptosporidium* (equivalent to 133 oocysts detected per test). Because the GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ test detects soluble antigen in fecal specimens in addition to cysts and oocysts, this LOD study represents an estimate of analytical sensitivity based on purified *Giardia* cysts and *Cryptosporidium* oocysts. Clinical specimens contain varying amounts of free antigen per *Giardia* cyst or *Cryptosporidium* oocyst.

REPRODUCIBILITY

A total of 22 fecal specimens were pre-characterized by commercially available predicate devices. The samples included 6 *Giardia*-positive specimens (3 mid-range positives), 6 *Cryptosporidium*-positive specimens (3 low positives), and 4 *Giardia/Cryptosporidium*-positive specimens (2 of which were low *Giardia*-positives and 2 of which were low *Cryptosporidium*-positives), and 6 specimens negative for both parasites. All specimens were coded to prevent their identification during testing. Testing was performed at 3 sites. The samples were tested, twice a day over a 5-day period by multiple technicians at each site using 2 different kit lots. A positive and negative control was run with each panel of the masked samples. The results from each laboratory were subsequently submitted to TECHLAB, Inc. and compared with in-house results. The results were consistent among the different locations, and exhibited a correlation of 100%. The positive specimens consistently tested positive and the negative specimens consistently tested negative at all sites using the GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ test.

CROSS REACTIVITY

The GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ test was evaluated for cross-reactivity with the bacterial and viral strains listed below. None of the strains were shown to cross-react with the GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™ test.

Aeromonas hydrophila
Bacillus subtilis

Bacillus cereus
Bacteroides fragilis

<i>Campylobacter coli</i>	<i>Campylobacter fetus</i>
<i>Campylobacter jejuni</i>	<i>Candida albicans</i>
<i>Clostridium bifermentans</i>	<i>Clostridium difficile</i> (strain 630)
<i>Enterococcus faecalis</i>	<i>Escherichia coli</i>
<i>Escherichia coli</i> O157:H7	<i>Escherichia coli</i> ETEC (enterotoxigenic)
<i>Escherichia coli</i> EPEC (enteropathogenic)	<i>Escherichia coli</i> EIEC (enteroinvasive)
<i>Klebsiella pneumoniae</i>	<i>Salmonella typhimurium</i>
<i>Shigella dysenteriae</i>	<i>Shigella flexneri</i>
<i>Shigella sonnei</i>	<i>Staphylococcus aureus</i>
<i>Staphylococcus aureus</i> (Cowan's)	<i>Staphylococcus epidermidis</i>
<i>Vibrio parahaemolyticus</i>	<i>Yersinia enterocolitica</i>
Human Adenovirus 1 and 3	Adenovirus, Type 2, 5, 40 and 41
Human Coxsackievirus B2, B3, and B4	Coxsackievirus B5
Human Coronavirus	Echovirus 11, 18, 33
Human Echovirus 9	Human paraechovirus 1 (Echovirus 22)
Enterovirus 68, 69	Human Enterovirus 70, 71
Human rotavirus	

Additionally, the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test was run on fecal specimens documented to be positive for other parasites by microscopy. No cross-reactivity was seen with the following organisms for either the *Giardia*-portion or the *Cryptosporidium*-portion of the test. Cross-reactivity to Astrovirus and Caliciviruses has not been established.

<i>Ascaris lumbricoides</i> eggs	<i>Blastocystis hominis</i>
<i>Chilomastix mesnili</i>	<i>Cyclospora cayetanensis</i>
<i>Dientamoeba fragilis</i>	<i>Diphyllobothrium latum</i> eggs
<i>Endolimax nana</i>	<i>Entamoeba coli</i>
<i>Entamoeba hartmanni</i>	<i>Entamoeba histolytica/E. dispar</i>
Hookworm eggs	<i>Iodamoeba bütschlii</i>
<i>Trichuris trichiura</i> eggs	

INTERFERING SUBSTANCES (U.S. Formulations)

The following substances had no effect on positive or negative test results analyzed at the concentrations indicated: Hog gastric mucin (3.5% w/v), Human blood (40% v/v), Barium sulfate (5% w/v), Imodium® (5% v/v), Kaopectate® (5% v/v), Pepto-Bismol® (5% v/v), Steric/Palmitic Acid (40% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v).

SPECIMEN PRESERVATIVE COMPATIBILITY

The *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* test was examined for compatibility with preserved fecal specimens during an in-house evaluation. No decrease in sensitivity or specificity was observed with specimens preserved in 10% buffered formalin, Sodium Acetate Formalin (SAF), Cary-Blair, or C&S when compared to fresh and frozen specimens.

PRECISION – INTRA-ASSAY

For the determination of intra-assay performance, 6 positive fecal specimens (two positive for *Giardia*, two positive for *Cryptosporidium*, two positive for both *Giardia* and *Cryptosporidium*) and six negative fecal specimens were analyzed. Each specimen was assayed on 5 cassettes. All positives remained positive and all negatives remained negative.

PRECISION – INTER-ASSAY

For the determination of inter-assay performance, 16 positive fecal specimens (six positive for *Giardia*, six positive for *Cryptosporidium*, and four positive for both *Giardia* and *Cryptosporidium*) and six negative fecal specimens were assayed twice a day over a four-day period using the *GIARDIA/CRYPTOSPORIDIUM QUIK CHEK™* kit. All positives remained positive and all negatives remained negative.

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Human Adenovirus 1 and 3
 Human Coxsackievirus B2, B3, and B4
 Human Coronavirus
 Human Echovirus 9
 Enterovirus 68, 69
 Human rotavirus

Adenovirus, Type 2, 5, 40 and 41
 Coxsackievirus B5
 Echovirus 11, 18, 33
 Human paraechovirus 1 (Echovirus 22)
 Human Enterovirus 70, 71

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Technical Support

Further information can be obtained from contacting TECHLAB® Technical Support::

US +1 800 TECHLAB
Phone (540) 953-1664
Fax (540) 953-1665
Email ts@techlab.com

© 2023 TECHLAB, Inc. All rights reserved.

GIARDIA/CRYPTOSPORIDIUM QUIK CHEK, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc.

All other trademarks referenced are trademarks of their respective owners.